**FRANCAIS**
**OptiMesh®**

#### DESCRIPTION

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**ENGLISH**
**OptiMesh®**

**DESCRIPTION**

- OptiMesh® is a sterile, conformable, three-dimensional mesh graft containment device knitted from polyester (polyethylene terephthalate (PET)) yarn. PET material is not bioabsorbable.
- OptiMesh is designed to contain and reinforce bone graft materials and/or bone graft substitutes, thereby limiting the migration of the graft materials into unintended areas.
- OptiMesh is offered in a range of shapes and sizes and is supplied attached to a disposable, delivery/holder tip for ease of handling and placement.
- OptiMesh is radiolucent, but may be printed with radiopaque reference markers. OptiMesh is cleaned after printing, packaged into double Tyvek/Mylar film pouches, placed in a carton and then terminally sterilized by gamma radiation.
- OptiMesh is provided in two pore sizes. The pore size variation permits vascular and bony ingrowth (1500 series) or provides restricted flow (500 series) of viscous bone substitute materials.
- OptiMesh is expandable/inflatable upon filling, requiring minimal exposure to position and deploy a relatively large device.
- Filling the OptiMesh device as directed, with use of the diverted tubes and granular graft, can generate a distraction/lifting force.
- OptiMesh is suitable for use with adjunct stabilizing fixation.

*Spineology Inc. expressly warrants that these devices are fabricated from the materials described. No other warranties, expressed or implied, are made.*

**INTENDED USE / INDICATIONS**

OptiMesh is a surgical mesh intended to contain and reinforce bone grafting materials and/or bone graft substitutes in bony voids or gaps of the skeletal system (i.e., the extremities, spine, and pelvis). These voids or gaps may be surgically created osseous defects, natural gaps, or defects partially surrounded by bone, defects created from traumatic injury, or a combination thereof. Adjunct fixation is recommended in cases involving segmental instability.

OptiMesh is indicated for use in orthopedic grafting procedures, including the following: Lumbar interbody fusion. The OptiMesh® 1500 System is used to contain, enhance the compaction of, and stabilize and reinforce bone graft material used in the intervertebral space of patients undergoing fusion to treat degenerative disc disease at one spinal level from L2-S1. For this use, OptiMesh 1500 is to be used in conjunction with stabilizing posterior fixation.

Thoracolumbar vertebral body fractures: OptiMesh® 1500 System with allograft is intended for vertebral augmentation, including fracture reduction, for the treatment of painful insufficiency fractures of the vertebral body due to osteoporosis, minor trauma or benign lesions. It is also intended for the treatment of painful traumatic vertebral fractures which are either a) clinically stable upon presentation, or b) are surgically stabilized with posterior fixation.

**CONTRAINDICATIONS**

- Infection, systemic or local, to the surgical site is a contraindication for placement of any surgical implant.
- Any medical condition that would preclude the patient from having surgery or would impede the benefit of implant surgery.
- Any patient unable or unwilling to adhere to physician's postoperative instructions or to participate in a postoperative management program (e.g., loss of musculature, neuromuscular compromise).
- Circulatory problems, such as thrombophlebitis, lymphedema, uncorrectable coagulopathy, or vascular deficiency at the implant site, when bone graft materials are used for filling.
- Patients requiring immediate subsequent radiation treatment for tumors, when bone graft materials are used for filling.

**PRECAUTIONS**

- As will any permanent implant, a perioperative antibiotic protocol is recommended.
- OptiMesh is part of a system. Spineology OptiMesh instruments are to be used as directed in the surgical technique manual for cavity preparation, proper mesh placement, and appropriate mesh filling.
- The surgeon should be trained and familiar with the OptiMesh System and with fluoroscopic image interpretation prior to surgery.
- OptiMesh is designed to work with granular allograft materials. Granules of non-demineralized corticocancellous bone provide structural properties for load-sharing support when contained by OptiMesh. Do not use demineralized bone pastes or putties alone to fill the OptiMesh device.
- OptiMesh is provided sterile. Product packaging should always be inspected for continuity and the components should be handled appropriately to ensure sterility. Do not use if packaging has been damaged or if the shelf life has been exceeded.
- OptiMesh is not intended to be trimmed or cut. The mesh can be damaged or torn through improper handling or usage. A torn mesh can lead to loss of graft containment.
- OptiMesh should never be reused.
- Resterilization of OptiMesh is not recommended under any circumstances.
- The OptiMesh Instrument Set is provided non-sterile and must be steam sterilized prior to use.
- Certain OptiMesh tools are provided as sterile single-use devices. Refer to the Surgical Technique Manual and package labels for additional information.

**GUIDELINES FOR USE**

- Use OptiMesh devices appropriate to the fill materials selected. The 500 series is intended for viscous bone substitutes (such as bone cement) and the 1500 series for bone graft (osteoconductive/osteoinductive) materials.
- Create or shape the cavity as required by the anatomy and pathology.

- Estimate the size of the cavity and select a suitable mesh.
- Place OptiMesh into prepared cavity.
- Fill OptiMesh using the filling tools appropriate for the mesh and fill material.
- Complete the procedure per the Surgical Technique Manual.

**CLEANING AND DECONTAMINATION**

**Cleaning Instructions**

- Cleaning and decontamination of surgical instruments is required before introduction into the sterile field.
- Following use, disassemble devices as instructed for cleaning.
- After use, preventing drying prior to cleaning will facilitate instrument cleaning.
- Soak in enzymatic detergent (mixed per manufacturer's recommendations) for one (1) minute or longer.
- Use a soft brush for manual cleaning and a soft bottle brush to clean the tube devices. Pay special attention to inner diameters and crevices during cleaning.
- Ultrasonic cleaning is acceptable, per established hospital methods.
- Rinse thoroughly under running water for one (1) minute or longer.

Surgical instruments must be handled with care. Improper handling may result in damage and may impair proper functioning of the device. Instruments which exhibit signs of damage or deterioration, including discoloration or corrosion, must be replaced. Ensure that all components of the system are available for use prior to surgery. Instruments must be sterilized before use and are to be cleaned and reesterilized immediately after use.

**STERILIZATION**

All instruments must be sterilized by the hospital before use and following any cleaning and decontamination procedures. Instrument trays are provided for storing and sterilizing the instruments. Sterilization can be performed on wrapped trays with the following cycle parameters:

Sterilization Method: Steam

Sterilizer Type/Cycle: Prevacuum

Temperature: 132 degrees C Minimum, 270 degrees F Minimum

Exposure Time: 4 Minutes Minimum

*Deviations from the recommended methods of cleaning and decontamination are not advised. It is the sole responsibility of the user to qualify such deviations.*

**FURTHER INFORMATION or PRODUCT COMPLAINTS:**

**Contact Spineology at:**

Spineology Inc.

St. Paul, MN 55128, USA

0120

**DEVICE RETRIEVAL:**

Contact Spineology to discuss Revision/Retrieval procedures.

**EC REP**
Quality First International
Bexhill-on-Sea, United Kingdom

**Catalog Description Information:**

OptiMesh pore size: 1500 micron pore size, for ingrowth applications: (1500E, 1500S, 1500DS); 500 micron pore size, for cement fill (500E).

The implant catalog numbers describe the group (first 3 digits), nomial filled body length (next 2 digits), and nominal neck length (final 2 digits).

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